mikrobiologie labor technik

# Chapman Agar

## PRINCIPIU

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
## PRINCIPLE

### Chapman Agar (Mannitol Salt Agar)

is used for the selective isolation and identification of Staphylococci in all types of samples. This medium is compliant with the formula described by Chapman, it contains a high level of Sodium Chloride (75 g/l) which inhibits the growth of most interfering flora (except a few *Bacillus*, *Corynebacterium* & *Streptococcus faecalis*). Acid productions resulting from the mannitol fermentation cause a colour change of the pH indicator, phenol red, from red to yellow. Pathogen Staphylococci, coagulase positive, grow as large colonies surrounded by a yellow halo. Staphylococci coagulase negative colonies are small and do not cause a change of colour.

## FORMULA

Components	g/l
Tryptone	5.00
Peptic digest of animal tissue	5.00
Meat extract	1.00
Mannitol	10.00
Sodium chloride	75.00
Phenol red	0.025
Agar	15.00
Final pH : 7.4 ± 0.2 at 25°C	

## METHOD

Suspend 111 grammes of powder in one litre of purified water. Bring to the boil under constant homogenisation until completely dissolved. Dispatch in tubes or flasks. Autoclave 15 minutes at 121°C.

## PROCEDURE

Liquefy the medium in boiling water then cool it to 45-50°C (according to standard

specifications) in a second water bath. Pour the medium into sterile Petri plates, let the plates set on an horizontal surface.

### Medical bacteriology:

Inoculate the surface of the plate with the specimen. Turn the plates upside-down before incubating them at 36 ± 2°C. Cultures are examined generally after 24 to 48 hours incubation.

### Harmonized method of pharmacopoeias:

Inoculate a plate of Chapman with the enriched sample. Turn the plate upside-down then incubate it at 32.5 ± 2.5°C for 18 to 72 hours.

### Water analyses: (According to a French standard)

Filter the appropriate volume of sample onto a membrane with 0.45µm pores. Transfer the membrane onto the Chapman plate. Take care not to capture air between the membrane and the agar. Turn the plate upside-down before incubating at 36 ± 2°C for 44 ± 4 hours.

## RESULTS

Pathogen Staphylococci grow as luxuriant colonies that produce their own pigment. These colonies are surrounded by a yellow halo due to mannitol's fermentation. Non pathogen Staphylococci grow mainly as small red colonies.

## REMARKS

The temperature of incubation of the agar is chosen according to procedure and current regulations. *Staphylococcus epidermidis* strains are able to ferment mannitol. Under those circumstances it is essential to carry out further investigation by look for example for the production of free coagulase and/or the A protein, clumping factor or even Dnase test.

## **PRECAUTIONS**

It is important not to dry the medium in an incubator since this might increase the sodium chloride concentration making the medium to selective.

## **BIBLIOGRAPHY**

1. Chapman G.H. 1945. The significance of sodium chloride in studies of staphylococci. J. Bact. 50:201-203.
2. Chapman G.H. 1948. An improved Stone medium for the isolation and testing of food poisoning staphylococci. Food Research, 13:100-105.
3. The United States Pharmacopeia.

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